

Anti-Biglycan Rabbit antibody



Catalog Number: 169087

Product name

Rabbit IgG

Anti-Biglycan Rabbit antibody

Applications

Specificity

WB, IP, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/500-1/2000

Biglycan Rabbit polyclonal antibody

IP: 1/50

Preparation

FC: 1/100

Antigen: A synthesized peptide derived from human Biglycan

Validations


Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Ig Type



Western blot analysis of Biglycan expression in C2C12 cell lysate.