

# HMGN2 antibody



Catalog Number: 111487

## Product name

HMGN2 antibody

## Immunogen

[Human HMGN2 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

HMGN2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, HepG2 cells, HL-60 cells, K-562 cells. Positive IHC detected in human prostate cancer tissue, human cervical cancer tissue. Observed molecular weight by Western-blot: 18-20 kDa

## Preparation

This antibody was obtained by immunization of HMGN2 recombinant protein (Accession Number: NM\_005517). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111487(HMGN2 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:111487(HMGN2 antibody) at dilution of 1:50 (under 10x lens)