

# Rb antibody

Catalog Number: 114487

---

## Product name

Rb antibody

## Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

## Antibody description

Rb Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells. Observed molecular weight by Western-blot: 110kd

## Preparation

This antibody was obtained by immunization of Rb recombinant protein (Accession Number: NM\_000321). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:114487(RB1 antibody) at dilution of 1:500