

Anti-CXCL3 antibody



Catalog Number: 176487

Product name

Polyclonal

Anti-CXCL3 antibody

Ig Type

Specificity

Rabbit IgG

Human

Applications

Antibody description

IHC-P

Rabbit polyclonal antibody to CXCL3

Dilutions

Preparation

IHC-P: 1:50-200

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human cxcl3. aa range:58-107

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1:: Immunohistochemical analysis of paraffin-embedded human-heart, antibody was diluted at 1:200