

# TUBGCP4 antibody



Catalog Number: 116487

## Product name

TUBGCP4 antibody

ELISA, WB, IHC

## Specificity

Human, Mouse, Rat; other species not tested.

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Antibody description

TUBGCP4 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue. Positive WB detected in MCF7 cells, HepG2 cells. Observed molecular weight by Western-blot: 76 kDa

## Validations




## Preparation

This antibody was obtained by immunization of TUBGCP4 recombinant protein (Accession Number: BC000966). Purification method: Antigen affinity purified.

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116487(TUBGCP4 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human placenta using Catalog No:116487(TUBGCP4 antibody) at dilution of 1:100 (under 10x lens)

## Applications