

PPAR alpha Polyclonal Antibody



Catalog Number: 160887

Product name

PPAR alpha Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PPAR alpha

Preparation

Antigen: Recombinant protein of human PPAR alpha

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IP


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IP 1:20 - 1:50

Validations



Western blot - PPAR alpha Polyclonal Antibody

Western blot analysis of extracts of U-937 cells, using PPAR alpha antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.