

# Piezo1 antibody



Catalog Number: 113887

---

## Product name

Piezo1 antibody

## Specificity

Human; other species not tested.

## Antibody description

Piezo1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, COLO 320 cells, HeLa cells, human brain tissue, HUVEC cells. Observed molecular weight by Western-blot: 230-280 kDa

## Preparation

This antibody was obtained by immunization of Piezo1 recombinant protein (Accession Number: BC008073). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113887(FAM38A antibody) at dilution of 1:400