

ALX3 antibody



Catalog Number: 107887

Product name

ALX3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ALX3 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of ALX3 recombinant protein (Accession Number: NM_006492). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107887(ALX3 Antibody) at dilution of 1:1000