

Dematin antibody



Catalog Number: 109887

Product name

Dematin antibody

Immunogen

[Human EPB49 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Dematin Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 35-45 kDa

Preparation

This antibody was obtained by immunization of Dematin recombinant protein (Accession Number: NM_001978). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109887(EPB49 Antibody) at dilution of 1:300