

NIP7 antibody

Catalog Number: 113187

Product name

NIP7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NIP7 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells, HepG2 cells. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 20,15kd

Preparation

This antibody was obtained by immunization of NIP7 recombinant protein (Accession Number: NM_016101). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113187(NIP7 antibody) at dilution of 1:500



IP Result of anti-NIP7 (IP:Catalog No:113187, 3ug; Detection:Catalog No:113187 1:500) with mouse liver tissue lysate 4000ug.