

Anti-Islet 1 Rabbit antibody



Catalog Number: 167187

Product name

Anti-Islet 1 Rabbit antibody

Specificity

Human, Mouse

Antibody description

Islet 1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Islet 1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of Islet 1 in Jurkat, K562 cell lysates using Islet 1 Rabbit pAb(1:1000 diluted). Predicted band size: 39KDa. Observed band size: 39KDa.