

Human PPIC/CYPC (Trx & His tag) recombinant protein



Catalog Number: 518187

General Information

Gene Name Synonym

Peptidyl-Prolyl Cis-Trans Isomerase C, PPIase C, Cyclophilin C, Rotamase C, PPIC, CYPC

Protein Construction

A DNA sequence encoding the Human PPIC/CYPC Lys31-Asp182 is expressed with a Trx, 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PPIC/CYPC (Trx & His tag) recombinant protein