

# CEP78 antibody

Catalog Number: 109187

---

**Product name**

CEP78 antibody

**Specificity**

Human; other species not tested.

**Antibody description**

CEP78 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells.

**Preparation**

This antibody was obtained by immunization of Peptide (Accession Number: NM\_032171).  
Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, IF

**Dilutions**

Recommended Dilution:

IF: 1:50-1:500

**Validations**

Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Catalog No:109187(CEP78 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)