

Human METAP1/KIAA0094 recombinant protein



Catalog Number: 519187

General Information

Gene Name Synonym

Methionine aminopeptidase 1; MAP 1; MetAP 1; Peptidase M 1; METAP1

Protein Construction

A DNA sequence encoding the Human METAP1/KIAA0094 Met1-Phe386 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 500mM NaCl, 10% glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human METAP1/KIAA0094 recombinant protein