

# Anti-Niemann Pick C2 Rabbit antibody



Catalog Number: 167587

## Product name

Rabbit IgG

Anti-Niemann Pick C2 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Niemann Pick C2 Rabbit monoclonal antibody

IHC: 1/500

## Preparation

ICC/IF: 1/200

Antigen: Recombinant protein of human Niemann Pick C2

FC: 1/20

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Niemann Pick C2 in HeLa, A549, HL-60 cell lysates using Niemann Pick C2 Rabbit mAb (1:1000 diluted). Predicted band size: 17kDa. Observed band size: 17kDa.