

# PODN antibody



Catalog Number: 113987

## Product name

PODN antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PODN Rabbit Polyclonal antibody. Positive WB detected in A375 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of PODN recombinant protein (Accession Number: XM\_011540677). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:113987(PODN Antibody) at dilution of 1:600