

# XIRP2 antibody

Catalog Number: 116987

## Product name

XIRP2 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

XIRP2 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle. Positive WB detected in COLO 320 cells, human heart tissue. Observed molecular weight by Western-blot: 260-300kd

## Preparation

This antibody was obtained by immunization of XIRP2 recombinant protein (Accession Number: BC022888). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116987(XIRP2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:116987(XIRP2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:116987(XIRP2 antibody) at dilution of 1:50

# XIRP2 antibody



Catalog Number: 116987

---

(under 40x lens)