

Anti-SIRT3 Rabbit antibody



Catalog Number: 168987

Product name

Anti-SIRT3 Rabbit antibody

Specificity

Human

Antibody description

SIRT3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human SIRT3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50

Validations



Western blot detection of SIRT3 in Jurkat, HeLa cell lysates using SIRT3 Rabbit mAb(1:500 diluted). Predicted band size: 44kDa. Observed band size: 28kDa.