

# EBP2 antibody



Catalog Number: 110287

## Product name

EBP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EBP2 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, A549 cells, HEK-293 cells, HeLa cells, NIH/3T3 cells, Raji cells. Observed molecular weight by Western-blot: 35-38 kDa

## Preparation

This antibody was obtained by immunization of EBP2 recombinant protein (Accession Number: NM\_006824). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110287(EBP2 Antibody) at dilution of 1:2000