

NSE antibody

Catalog Number: 113287

Product name

NSE antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NSE Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive WB detected in SH-SY5Y cells, HEK-293 cells, HepG2 cells, mouse brain tissue. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001975). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113287(ENO2 antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:113287(ENO2 Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:113287(ENO2 Antibody) at dilution of 1:200 (under 40x lens)