

# PSAP antibody



Catalog Number: 114287

## Product name

PSAP antibody

Recommended Dilution:

WB: 1:200-1:1000

## Specificity

Human, Mouse; other species not tested.

IHC: 1:20-1:200

## Antibody description

PSAP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, MCF7 cells, mouse brain tissue. Positive IHC detected in human lung cancer tissue, human brain tissue, human kidney tissue, human liver tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of PSAP recombinant protein (Accession Number: NM\_002778). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


ELISA, IHC, WB

## Dilutions


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114287(PSAP antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114287(PSAP antibody) at dilution of 1:200 (under 10x lens)



# PSAP antibody



Catalog Number: 114287

---

Immunohistochemical of paraffin-embedded  
human lung cancer using Catalog

No:114287(PSAP antibody) at dilution of 1:200  
(under 40x lens)