

# SIGLEC5 antibody



Catalog Number: 115287

## Product name

SIGLEC5 antibody

## Specificity

Human; other species not tested.

## Antibody description

SIGLEC5 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, K-562 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of SIGLEC5 recombinant protein (Accession Number: XM\_047446914). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:115287(SIGLEC5 antibody) at dilution of 1:500