

# Mouse anti mCherry-Tag mAb



Catalog Number: 166287

## Product name

Mouse anti mCherry-Tag mAb

## Antibody description

Monoclonal antibody to Mouse anti mCherry-Tag

## Preparation

Antigen: Recombinant protein of human mCherry-tag

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications


WB, IF

## Dilutions

WB 1:2000 - 1:5000

IF 1:50 - 1:200

## Validations



Western blot - Mouse anti mCherry-Tag mAb

Western blot analysis of over-expressed mCherry-tagged protein in 293T cells using Mouse anti mCherry-Tag mAb at different dilution.