

Anti-CD95/APO-1/TNFRSF6 antibody



Catalog Number: 106287

Product name

Anti-CD95/APO-1/TNFRSF6 antibody

Immunogen

[Rat CD95/APO-1/TNFRSF6 \(His Tag\) recombinant protein](#)

Specificity

Rat CD95 / APO-1 / TNFRSF6 / FAS

Antibody description

Rabbit monoclonal to CD95/APO-1/TNFRSF6

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat CD95 / APO-1 / TNFRSF6 / FAS (rR CD95 / APO-1 / TNFRSF6 / FAS; NP_631933.2; Met1-Lys170).

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM: 0.5-2 μ g/Test

Validations



CD95 / APO-1 / TNFRSF6 / FAS Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Rat FAS(CD95) expression on SD rat splenocytes. Cells were stained with purified anti-Rat FAS(CD95), then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.