

Human APOC2/APC2 (His tag) recombinant protein



Catalog Number: 518287

General Information

Gene Name Synonym

Apolipoprotein C-II, Apolipoprotein C2, APC2 and APOC2

Protein Construction

A DNA sequence encoding the Human APOC2/APC2 Met1-Glu82 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mMPB, 150mM NaCL, 30%glycerol, pH7.2.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human APOC2/APC2 (His tag) recombinant protein