

Anti-CD22 antibody



Catalog Number: 104687

Product name

Anti-CD22 antibody

Immunogen

[Rat CD22 \(His Tag\) recombinant protein](#)

Specificity

Rat CD22

Antibody description

Rabbit polyclonal to CD22

Preparation

Produced in rabbits immunized with purified, recombinant Rat CD22 (rR CD22; NP_001100973.1; Met1-Arg692). CD22 specific IgG was purified by Rat CD22 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Rat CD22.

IHC-P: 0.1-2 µg/mL

Validations



CD22 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of rat CD22 in rat spleen with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).