

Anti-UTF1 antibody

Catalog Number: 176687

Product name

Anti-UTF1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to UTF1

Preparation

This antigen of this antibody was peptide

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:1,000-1:2,000

IHC-P: 1:200

Validations




Fig1: Western blot analysis on D3 cell lysates (1) and human brain tissue (2) using anti-UTF1 rabbit polyclonal antibodies.




Fig2: Immunohistochemical analysis of paraffin-embedded human testis tissue using anti- UTF1 rabbit polyclonal antibodies.