

# TTC7A antibody

Catalog Number: 116387

## Product name

TTC7A antibody

WB: 1:200-1:1000

IP: 1:200-1:1000

## Specificity

Human, Mouse; other species not tested.

IHC: 1:20-1:200

## Antibody description

TTC7A Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, mouse colon tissue, mouse small intestine tissue, PC-3 cells, Transfected HEK-293 cells. Positive IHC detected in human colon tissue. Observed molecular weight by Western-blot: 96 kDa

## Preparation

This antibody was obtained by immunization of TTC7A recombinant protein (Accession Number: NM\_001288955). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IP


## Dilutions

Recommended Dilution:


## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116387(TTC7A Antibody) at dilution of 1:300



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116387(TTC7A Antibody) at dilution of 1:500




# TTC7A antibody

Catalog Number: 116387


---

IP Result of anti-TTC7A (IP:Catalog No:116387, 4ug; Detection:Catalog No:116387 1:300) with K-562 cells lysate 4000ug.



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog

No:116387(TTC7A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:116387(TTC7A Antibody) at dilution of 1:50 (under 40x lens)