

# Anti-SIRT1 (2G10) Mouse antibody



Catalog Number: 168387

## Product name

Anti-SIRT1 (2G10) Mouse antibody

## Specificity

Human

## Antibody description

SIRT1 (2G10) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human SIRT1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB, IP


## Dilutions

WB: 1/1000

## Validations



Western blot detection of SIRT1 in HeLa, Jurkat, Ramos, A549 and CEM cell lysates using SIRT1 mouse mAb (1:1000 diluted). Predicted band size: 120KDa. Observed



band size: 120KDa.

Immunoprecipitation analysis of HeLa cell lysates using SIRT1 mouse mAb.