

# FTSJ2 antibody

Catalog Number: 110787

## Product name

FTSJ2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FTSJ2 Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in human skeletal muscle tissue, human placenta tissue. Observed molecular weight by Western-blot: 29 kDa

## Preparation

This antibody was obtained by immunization of FTSJ2 recombinant protein (Accession Number: NM\_013393). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110787(FTSJ2 Antibody) at dilution of 1:600



IP Result of anti-FTSJ2 (IP:Catalog No:110787, 3ug; Detection:Catalog No:110787 1:500) with mouse skeletal muscle tissue lysate 4000ug.