

MTMR4 antibody



Catalog Number: 112787

Product name

MTMR4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MTMR4 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, HepG2 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 140-160kd

Preparation

This antibody was obtained by immunization of MTMR4 recombinant protein (Accession Number: NM_004687). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:112787(MTMR4 antibody) at dilution of 1:500