

VPS39 antibody



Catalog Number: 116787

Product name

VPS39 antibody

Specificity

Human; other species not tested.

Antibody description

VPS39 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 102 kDa

Preparation

This antibody was obtained by immunization of VPS39 recombinant protein (Accession Number: BC015817). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-VPS39 (IP:Catalog No:116787, 4ug; Detection:Catalog No:116787 1:300) with mouse brain tissue lysate 4000ug.