

# Cytohesin 1 antibody



Catalog Number: 109787

## Product name

Cytohesin 1 antibody

## Immunogen

[Human CYTH1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Cytohesin 1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of Cytohesin 1 recombinant protein (Accession Number: NM\_004762). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109787(CYTH1 antibody) at dilution of 1:500