

DNAJB6 antibody

Catalog Number: 110097

Product name

DNAJB6 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Validations

Antibody description

DNAJB6 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue, HEK-293 cells, Y79 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 27kd



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110097(DNAJB6 antibody) at dilution of 1:100 (under 10x lens)

Preparation

This antibody was obtained by immunization of DNAJB6 recombinant protein (Accession Number: NM_005494). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, IF, WB, IP

Dilutions

Recommended Dilution:



WB: 1:200-1:2000



Immunofluorescent analysis of HepG2 cells, using DNAJB6 antibody Catalog No:110097 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

DNAJB6 antibody

Catalog Number: 110097



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110097(DNAJB6 Antibody) at dilution of 1:600

IP Result of anti-DNAJB6 (IP:Catalog No:110097, 4ug; Detection:Catalog No:110097 1:800) with mouse brain tissue lysate 4000ug.