

Anti-ACT11 antibody



Catalog Number: 175097

Product name

Anti-ACT11 antibody

Specificity

Arabidopsis thaliana, Oryza sativa, Nicotiana tabacum

Antibody description

Mouse monoclonal antibody to ACT11

Preparation

This antigen of this antibody was arabidopsis act11 full-length protein.

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, 50%Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store the product at -20°C. Stable for one year from the date of shipment.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

WB

Dilutions

WB:1:1,000-1:5,000

Validations




Fig1.; Western blot analysis of Actin on Nicotiana tabacum (1) and Oryza sativa (2) lysates using anti-Actin antibody.