

Anti-NOX2/gp91phox Rabbit antibody



Catalog Number: 167097

Product name

Polyclonal

Anti-NOX2/gp91phox Rabbit antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat

Applications

Antibody description

WB

NOX2/gp91phox Rabbit polyclonal antibody

Dilutions

Preparation

WB: 1/1000-1/2000

Antigen: A synthesized peptide derived from human CYBB

Validations


Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality



Western blot analysis of NOX2/CYBB/gp91phox expression in MCF-7 cell lysate.