

Anti-FH/Fumarase (9B1) Mouse antibody



Catalog Number: 168097

Product name

Anti-FH/Fumarase (9B1) Mouse antibody

WB: 1/3000

IF: 1/200

Specificity

Human, Mouse, Rat

IHC: 1/50-300

Antibody description

FH/Fumarase (9B1) Mouse monoclonal antibody

Preparation

Antigen: Synthetic Peptide of FH

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Mouse IgG1

Applications


WB, IHC-P, ICC/IF

Dilutions

Validations



Western blot detection of Fumarase in MCF7, Jurkat and HeLa cell lysates using Fumarase mouse mAb (1:3000 diluted). Predicted band size: 50 kDa. Observed band size: 50 kDa.



Western blot detection of Fumarase in K562, Ramos, Jurkat, HeLa cell lysates using Fumarase (9B1) Mouse mAb (1:1000 diluted). Predicted band size: 50 kDa. Observed band size: 50 kDa.