

# Anti-FH/Fumarase (9B1) Mouse antibody



Catalog Number: 168097

## Product name

Anti-FH/Fumarase (9B1) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FH/Fumarase (9B1) Mouse monoclonal antibody

## Preparation

Antigen: Synthetic Peptide of FH

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/3000


IF: 1/200

IHC: 1/50-300

## Validations



Western blot detection of Fumarase in MCF7, Jurkat and HeLa cell lysates using Fumarase mouse mAb (1:3000 diluted). Predicted band size: 50 kDa. Observed band size: 50 kDa.



Western blot detection of Fumarase in K562, Ramos, Jurkat, HeLa cell lysates using Fumarase (9B1) Mouse mAb (1:1000 diluted). Predicted band size: 50 kDa. Observed band size: 50 kDa.