

Mouse beta-Catenin/CTNNB1 (His & GST Tag) recombinant protein



Catalog Number: 500497

General Information

Gene Name Synonym

Beta-catenin

Protein Construction

A DNA sequence encoding the mouse CTNNB1 (Q02248) (Met1-Leu781) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Mouse

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant mouse CTNNB1/GST chimera

consists of 1018 amino acids and has a calculated molecular mass of 113 kDa. The recombinant protein migrates as an approximately 115 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 8.0, 10% gly

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse beta-Catenin / CTNNB1 Protein (His & GST Tag) SDS-PAGE