

SOHLH1 antibody



Catalog Number: 115497

Product name

SOHLH1 antibody

Specificity

Human; other species not tested.

Antibody description

SOHLH1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of SOHLH1 recombinant protein (Accession Number: NM_001012415). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115497(SOHLH1 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115497(SOHLH1 Antibody) at dilution of 1:100 (under 40x lens).