

Pan Phospho-Tyrosine mAb



Catalog Number: 166497

Product name

Pan Phospho-Tyrosine mAb

Specificity

ALL

Antibody description

Monoclonal antibody to Pan Phospho-Tyrosine

Preparation

Antigen: Phosphotyramine-KLH

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB, IHC, IP


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IP 1:50 - 1:100

Validations



Western blot - Pan Phospho-Tyrosine mAb

Western blot analysis of extracts of THP-1 cells, using Pan Phospho-Tyrosine antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.