

Anti-PDE1B Rabbit antibody



Catalog Number: 167497

Product name

Anti-PDE1B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PDE1B Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human PDE1B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/20

Validations



Western blot detection of PDE1B in Rat Brain lysates using PDE1B Rabbit pAb(1:1000 diluted). Predicted band size:61kDa.Observed band size:61kDa.