

COX7A1 antibody



Catalog Number: 109497

Product name

COX7A1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

COX7A1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse heart tissue, human brain tissue, human kidney tissue. Observed molecular weight by Western-blot: 7-9kd

Preparation

This antibody was obtained by immunization of COX7A1 recombinant protein (Accession Number: NM_001864). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109497(COX7A1 antibody) at dilution of



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109497(COX7A1 antibody) at dilution of 1:100 (under 10x lens)