

POMGNT2 Polyclonal Antibody



Catalog Number: 164897

Product name

POMGNT2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to POMGNT2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 301-580 of human POMGNT2 (NP_116195.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:20 - 1:50

Validations



Immunofluorescence - POMGNT2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using POMGNT2 antibody . Blue: DAPI for nuclear staining.