

Anti-CD2BP2 antibody

Catalog Number: 175897

Product name

Anti-CD2BP2 antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to CD2BP2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cd2 tail-binding. aa range:101-150

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using CD2BP2 Polyclonal Antibody diluted at 1: 1000