

# TCEA1 antibody

Catalog Number: 115897

## Product name

TCEA1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TCEA1 Rabbit Polyclonal antibody. Positive WB detected in Sp2/0 cells, HeLa cells, SGC-7901 cells. Positive IP detected in HeLa cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 40 kDa

## Validations




## Preparation

This antibody was obtained by immunization of TCEA1 recombinant protein (Accession Number: NM\_006756). Purification method: Antigen affinity purified.

Sp2/0 cells were subjected to SDS PAGE followed by western blot with Catalog No:115897(SII-TFIIS, TCEA1 antibody) at dilution of 1:800

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:115897(SII-TFIIS, TCEA1 antibody) at dilution of 1:400 (under 10x lens)

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000




Immunohistochemical of paraffin-embedded

# TCEA1 antibody

Catalog Number: 115897

human brain using Catalog No:115897(SII-TFIIS,TCEA1 antibody) at dilution of 1:400 (under

IP Result of anti-SII-TFIIS,TCEA1 (IP:Catalog No:115897, 4ug; Detection:Catalog No:115897 1:500) with HeLa cells lysate 520ug.



40x lens)