

Anti-AMPK beta 1 Rabbit antibody



Catalog Number: 166897

Product name

Anti-AMPK beta 1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

AMPK beta 1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human AMPK beta 1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/20-1/500

ICC/IF: 1/100

FC: 1/20

IP: 1/20

Validations



Western blot detection of AMPK β 1 in HeLa,PC12,3T3,COS7 cell lysates using AMPK β 1 Rabbit pAb(1:1000 diluted).Predicted band size:30KDa.Observed band size:38KDa.