

YWHAB antibody



Catalog Number: 116897

Product name

YWHAB antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

YWHAB Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue.

Preparation

This antibody was obtained by immunization of YWHAB recombinant protein (Accession Number: BC010352). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116897(YWHAB antibody) at dilution of 1:50 (under 10x lens)