

Anti-LYRIC/AEG1 Rabbit antibody



Catalog Number: 167897

Product name

Rabbit IgG

Anti-LYRIC/AEG1 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/2000-1/10000

LYRIC/AEG1 Rabbit monoclonal antibody

IHC: 1/50

Preparation

ICC/IF: 1/100

Antigen: A synthetic peptide of human LYRIC/AEG1

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of LYRIC/AEG1 in CHO-K1, Ramos using LYRIC/AEG1 Rabbit mAb(1:1000 diluted)