

Advillin antibody



Catalog Number: 107897

Product name

Advillin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Advillin Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 92kd

Preparation

This antibody was obtained by immunization of Advillin recombinant protein (Accession Number: BC111730). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:107897(AVIL antibody) at dilution of 1:300